

I Care Yal 240ml

Brand: ESOFORM
Product Code: 0018
Availability: In Stock
Weight: 0.00kg

Description

NEW PRODUCT

I Care Yal

We inform you of the new MULTI-PURPOSE CONTACT LENS CLEANING SOLUTION **I Care Yal** (disinfection - maintenance - cleaning - rinsing - hydration of all soft contact lenses - protein removal),

launching with the following characteristics:

-Contains Sodium Hyaluronate

-Works within 5-7 minutes

-Recommended for sensitive eyes.

-Moisturizes the lenses, immediately after insertion and throughout the day, prolonging their comfortable use and reducing dryness and irritation of the eyes. Its PH and osmoticity have been specially studied to be completely biocompatible with ocular tissues and organic fluids.

-240 ml bottle solution content , expiring after 3 years of its production.

-Must be used within 3 months, after the first day of use.

-Free samples offered, upon request. (contact : - e-mail: theeyewearcompany@gmail.com -t: +30 2310 951 505)