

# BioDrop MD®

**Brand:** PILOSET

**Product Code:** 0031

**Availability:** In Stock

**Weight:** 0.02kg

## Short Description

The new **BioDrop MD®** comfort eye drop maintains eyes moist and relaxed. **Hyaluronan**, the natural moisturizer, aids symptoms of the **dry eye** disease, protects epithelium of the cornea and expedites healing of corneal wounds.

## Description

The new **BioDrop MD®** comfort eye drop maintains eyes moist and relaxed. **Hyaluronan**, the natural moisturizer, aids symptoms of the **dry eye** disease, protects epithelium of the cornea and expedites healing of corneal wounds.

To ensure eyes with a natural feeling, the chemical and physical properties of **BioDrop MD®** are adjusted to resemble the healthy tear film.

**BioDrop MD®** contains **neither preservatives nor ingredients** of animal or plant origin.

**BioDrop MD®** is suitable for frequent and prolonged use. It is packed into a handy multidose dispensing system (10 ml) designed for **sterile and preservative-free** ophthalmic products.

Developed by **Finnsusp** in co-operation with the Finnish Allergy and Asthma Federation. As with all other **Piiloset** products, **BioDrop MD®** is not animal tested.

## PIILOSET BioDrop MD®

Preservative-free multidose comfort eye drop

### NATURAL MOISTURE FOR ALL DRY EYES

The new BioDrop MD® comfort eye drop maintains eyes moist and relaxed. Hyaluronan, the natural moisturizer, aids symptoms of the dry eye disease, protects epithelium of the cornea and expedites healing of corneal wounds<sup>1,2</sup>. To ensure eyes with a natural feeling, the chemical and physical properties of BioDrop MD® are adjusted to resemble the healthy tear film.



### GOOD FOR SENSITIVE EYES

BioDrop MD® is one of the eye-friendly PiiLOSET Allergy Label products. Examined by specialists of the Finnish Allergy and Asthma Federation, BioDrop MD® is confirmed free from generally irritating ingredients. It contains neither preservatives nor ingredients of animal or plant origin.

Hyaluronan has been produced for BioDrop MD® by biotechnological manufacturing methods and extensively assessed for safety and chemical purity.



Finnsusp Ltd | Paaskyalliontie 5, 21420 Lieto, Finland | Tel +358 2 489 110 | info@piiloseit.fi | www.piiloseit.fi

### MORE PLEASURE FROM CONTACT LENSES

To take full advantage of the moisturizing PiiLOSET care products, we recommend using moisturizing BioTwin® contact lens care solution together with BioDrop MD® comfort drops while wearing contact lenses.

BioDrop MD® is suitable for frequent and prolonged use. It is packed into a handy squeezable multidose dispensing system designed for sterile and preservative-free ophthalmic products. The system is calibrated exclusively for the properties of BioDrop MD®.

### AVAILABILITY

10 ml multidose dispenser  
Shelf life 2 years

### DESIGNED FOR PRESERVATIVE FREE FORMULATIONS

#### Tip seal & filter

- THE TIP SEAL MECHANISM CLOSSES THE ORIFICE IMMEDIATELY WHEN RELEASING THE PRESSURE, THUS PREVENTING CONTAMINATION THROUGH THE TIP
- SUITABLE FOR OXIDIZING FORMULATIONS DUE TO METAL FREE FLUID PATH



- THE AIR REQUIRED TO COMPENSATE THE VOLUME LOSS IN THE BOTTLE AFTER THE ACTUATION IS FILTERED THROUGH A PTFE MULTILAYER MEMBRANE



<sup>1</sup>Argona, P. et al. Long-term treatment with sodium hyaluronate containing artificial tears reduces ocular surface damage in patients with dry eye. *Br J Ophthalmol* 2002; 86: 818-824.  
<sup>2</sup>Lee, JH et al. Efficacy of sodium hyaluronate and carbomethylcellulose in treating mild to moderate dry eye disease. *Cornea* 2011; 30: 1759-1763.